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- Cartridge(50 mL) x 4 ea - Cartridge(50 mL) x 10 ea

Please Note

- Working time includes mixing time.
- There could be about 15 seconds difference(faster or later) of working time per one Degree C. temperature change.
- Latex glove may disturb hardening of the impression materials, please use well-washed hands or vinyl glove
- Eugenol or hemostatic can also disturb hardening. If Hydrogen peroxide is used, wash out with warm water and eliminate bubbles.
- Keep away from the cloth, you can not remove it once it touches cloth.
- Pay attention not to change any cap or measuring spoons of each materials.
- Do not use for the patients who are allergic to silicon.
- There products are for dental clinic purpose only.
- Only professionals can handle these products.

Storage

- If the room temperature is over 25°C, store in cool place.
- Validity is 2 years from the date of manufacture.
- Defective products would be replaced free of charge.

<Sterilization>

- After remove the impression materials from the mouth, make sure to wash it. Keep 30 Min. for sterilization in the 2% glutaraldehyde(acidic), 0.5~1.5% Sodium hypodivone or 0.1~1.5% Providone-iodin.

<Plaster Pouring>

- Pour the plaster at least 30 Minutes after removal of the impression materials from the mouth for the best plaster model.

*Recovery time for the impression material : about 30 minutes.

* H2 Gas development time : 5-15minutes .

Made in Korea

Manufacturer	
Authorized Representative In the European Community	
Batch Code	
Date of Manufacture	
Use By	
Non Sterile	
Do not use if package is damaged	
Consult Instructions for use	
Temperature Limitation	
Keep dry	
Keep away from sunlight	
Caution	
Do not re-use	

B&E KOREA CO., LTD.
B&E KOREA | 취비엔이코리아
www.bnekorea.com

Manufacturer
B&E KOREA Co., Ltd.
#807, #802, B Dong, 60, Haan-ro,
Gwangmyeong-si, Gyeonggi-do, Korea
Tel. 82-70-7425-2818 Fax. 82-2-6112-2816

EC Representative
BIOGERI
Am Pfelshof 12, 22393 Hamburg, Germany
Tel. +49-40-4143-102150 Fax. +49-40-4143-102159

B&E KOREA CO., LTD.

Chemi SiL
Instruction for Use

Chemi Sil is hydrophilic
impression material based
on A-silicones

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Chemi SiL

Chemi Sil is hydrophilic impression material based on A-silicones

Characteristics

- Easy to mix
- High tensile strength
- Hydrophilic
- Very fine detail accuracy
- Homogeneous bubble free

Indication

- * Crown and bridge impressions, inlay and onlay impressions
- * Impression of the occlusal surface (Chemi SiL Bite)

* The matrix for Light Curing

Usage

- Insert the cartridge into the dispenser.
- Squeeze little bit until base and catalyst drain together without blocking.
- Insert the mixing tip into the cartridge and turn it 90 degree to the clock wise direction.
- Pull the trigger slowly for good mixing of the materials.

Working Time/setting Time

Products	Working time (23°C)	Setting time(Oral) (32°C)	Hardness (Shore A)
Chemi SiL LB/LBS(Light Body)	1 : 30	3 : 30	52 / 48
Chemi SiL MB/HB(Heavy Body)	1 : 30	3 : 30	70 / 60
Chemi SiL Bite-Blu	0 : 20	1 : 30	91
Sil Clear	0 : 30	2 : 00	60