Telio® Inlay / Onlay

[en] Instructions for Use

Light-curing temporary restorative material

Rx ONLY

CE 0123

Date information prepared:
2023-05-12/Rev.0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

- Temporary restoration of anterior and posterior teeth
- Temporary sealing of implant screw access openings in anterior and posterior teeth

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants

Intended users / Special training

- Dentists
- No special training required.

Use

For dental use only.

Description

Telio® Inlay and Telio® Onlay are light-curing materials for temporary restorations in anterior and posterior teeth and for temporary sealing of implant screw access openings in anterior and posterior teeth. These single-component materials for temporary restorations do not require the application of a temporary cement.

Both materials are available in a transparent (Trans.) and universal (Uni) shade and are supplied in syringes and Cavifils.

Indications

Missing tooth structure in anterior and posterior teeth

Areas of application (recommendation Telio Inlay):

- Temporization of deep inlay preparations with parallel walls
- Sealing of implant screw access holes

Areas of application (recommendation Telio Onlay):

- Temporization of large, shallow cavities (onlay preparations). Please note that after curing, Telio Onlay is less elastic than Telio Inlay, which may complicate removal if undercuts are present.
- Sealing of implant screw access holes

Contraindications

Telio Inlay and Telio Onlay must not be used:

- If the patient is known to be allergic to any of the ingredients of Telio Inlay or Telio Onlay

Limitations of use

- Not intended for the fabrication of temporary crowns and bridges
- Not intended for the fabrication of long-term provisional restorations (worn longer than 6 weeks)
-  Do not reuse the Cavifil. Cavifils must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).
- Do not use the syringe intraorally.

Processing restrictions

Do not use oxidizing disinfectants to disinfect syringes and Cavifils.

Side effects

After longer periods of service, the material may demonstrate discolouration. However, this does not compromise the performance of the material. In rare cases, Telio Inlay and Telio Onlay may cause a sensitizing reaction in patients with a hypersensitivity to any of the materials' ingredients. In these cases, the materials should no longer be used.

Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of Telio Inlay and Telio Onlay. Therefore, Telio Inlay and Telio Onlay should not be applied in combination with zinc oxide-eugenol cements.

Hemostatic substances may inhibit polymerization and/or lead to discolouration. Therefore, the directions for use of these substances must be observed.

Peroxide-based bleaching agents and disinfectants may inhibit the polymerization of Telio Inlay/Onlay. Such materials must not be used prior to the application of Telio Inlay/Onlay.

Clinical benefit

- Reconstruction of chewing function
- Temporary sealing of dentin lesions
- Restoration of esthetics

Composition

- Telio Inlay
Silicon dioxide, copolymer, polyester dimethacrylate, AAEMA, D3MA, potassium hexafluorotitanate
- Telio Onlay
Silicon dioxide, copolymer, polyester dimethacrylate, ETGMA, AAEMA, D3MA, potassium hexafluorotitanate

2 Application

2.1 Temporary restoration of anterior and posterior teeth

2.1.1 Cleaning and preparation

Clean cavity with water spray. Establish an adequate isolation of the prepared tooth; use a matrix if necessary.

2.1.2 Application

Place a sufficient quantity of Telio Inlay or Telio Onlay into the cavity using a spatula or another suitable instrument. If you use Cavifils, the material can be applied directly into the cavity by means of an application instrument. Spread out the material towards the preparation margin and remove any excess.

2.1.3 Light-curing of the temporary restoration

Increments of up to 4 mm can be cured in 10 seconds per surface with a standard curing light (at approx. 650 mW/cm²).

2.1.4 Finishing of the temporary restoration

- Remove matrix if required.
- Check occlusion and adjust if required.

2.1.5 Removal of the temporary restoration

- Insert a suitable instrument (e.g. probe/scaler) into the temporary restoration. Remove Telio Inlay or Telio Onlay from the cavity.
- Subsequently, clean the cavity with an oil- and fluoride-free cleaning paste.
- Try in the final restoration and seat it as usual.

2.2 Temporary sealing of implant screw access openings in anterior and posterior teeth

2.2.1 Cleaning

Clean screw channel with water spray. Establish an adequate isolation of the screw channel.

2.2.2 In order to ensure that the abutment screw is not

contaminated by temporary restorative material, apply teflon tape or a cotton pad in the screw channel.

Place a sufficient quantity of Telio Inlay or Telio Onlay into the screw channel using a spatula or another suitable instrument. If you use Cavifils, the material can be applied directly into the screw channel by means of an application instrument. Spread out the material towards the edge and remove any excess.

2.2.3 Light-curing of the temporary restoration

Increments of up to 4 mm can be cured in 10 seconds per surface with a standard curing light (at approx. 650 mW/cm²).

- 2.2.4 Finishing of the temporary restoration
Check occlusion and adjust if required.
- 2.2.5 Removal of the temporary restoration
 - Insert a suitable instrument (e.g. probe/scaler) into the temporary restoration. Remove Telio Inlay or Telio Onlay from the screw channel.
 - Rinse the screw channel with a stream of water and dry it with oil-free air.

Notes

- Telio Inlay / Telio Onlay remains well in place in standard preparations. By placing a wooden wedge prior to the application of the material, tighter proximal contacts and thus improved retention can be achieved.
- To make adjustments by grinding or remove excess material, use suitable polishers or tungsten carbide finishers. Excess material may also be removed with a scalpel or a scaler.
- Telio Inlay / Telio Onlay can be contoured more easily by wetting the instrument with an unfilled bonding agent.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available on the website www.ivoclar.com
- Explanation of the symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL014JK

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Avoid contact of uncured material with the skin, mucous membrane and eyes. Uncured material may cause skin irritations and may lead to sensitization. Commercially available medical gloves do not offer protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Loss of the filling may cause accidental ingestion of material and may make additional dental treatments necessary.
- Fracture of the filling may make additional dental treatments necessary.

4 Shelf life and storage

- Do not use Telio Inlay / Telio Onlay after the indicated expiry date.
- Expiry date: see information on syringes, Cavifils and packages. Do not use the product after the indicated expiry date.
- Close syringes/Cavifils immediately after use. Exposure to light leads to premature polymerization.
- Storage temperature 2–28 °C (36–82 °F)
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

For dental use only.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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