



TRIBEST DENTAL PRODUCTS CO., LTD

No.98 Zhongling Road, Sanmao Street, Yangzhong City, Jiangsu, China



Tel: +86-511-88091812

Fax: +86-511-88222305

Website: www.tribestdental.com

Email: kevin@tribest.cn

TECHNICAL DATA SHEET

Product Name	Self-Sealing Sterilization Pouch						
							
Composition	<p>blue tinted film sealed to paper with double-sides coated adhesive tape with two types of indicators: Steam: the indicator turns from blue to cocoa/brown when steam processed EO: the indicator turns from pink to blue when EO processed</p>						
Description	<p>Utilization: provide health care workers with an effective method to enclose devices intended for sterilization in steam autoclaves or via Ethylene Oxide (EO). The pouch's process indicators are designed to indicate to the user that the pouch has undergone either a steam or EO sterilization process.</p> <p>The pouches are made with high quality porous medical grade paper and a laminated PET/PP thick blue tinted film which ensure proper sterilant penetration.</p> <p>Pouches and packaging are free of lead/latex.</p> <p>Non-sterile</p>						
Size and Packaging (customizable)							
Single side indicator	11-50001	11-50002	11-50003	11-50004	11-50005	11-50006	
	57x130mm	70x260mm	90x165mm	90x260mm	135x260mm	135x290mm	
Double side indicator	11-50007	11-50008	11-50009	11-50010	11-50001	11-50001	
	57x130mm	70x260mm	90x165mm	90x260mm	135x260mm	135x290mm	
Packaging	200pcs/box 60boxes/ctn	200pcs/box 25boxes/ctn	200pcs/box 10boxes/ctn	200pcs/box 20boxes/ctn	200pcs/box 10boxes/ctn	200pcs/box 10boxes/ctn	
Comment	<p>We are ISO 13485 and CE certificated.</p> <p>Single use product. Do not reuse.</p> <p>Meets or exceeds EN868 standard.</p> <p>Shelf life: The pouches should be used within 3 years of manufacture. The sterilization pouch maintains the enclosed devices sterile up to minimum 6 months after sterilization.</p> <p>Storage: Keep away from light, heat and humidity. Store at temperature between 15 – 35°C with a relative humidity of <70% .</p>						
Manufacturer	<p>TRIBEST DENTAL PRODUCTS CO., LTD No.98 Zhongling Road, Sanmao Street, Yangzhong City, Jiangsu, China</p>						